

BRAUN

Nasal aspirator 1


Quick-start
guide

Guía de inicio
rápido





1 Insert Batteries
Inserte las pilas

2 x AA
(LR6 1.5V)





2 Select tip & press in place
Seleccione la punta y presione en
su lugar



High
Alto
Low
Bajo






Do not overfill
reservoir
No sobrellene
el depósito





Cleaning Instructions:
Twist to lock/unlock reservoir
Instrucciones de Limpieza:
Gire para asegurar/liberar el depósito



US For additional information visit:
www.braunnasalaspirator.com
SP Para obtener más información visite:
www.braunnasalaspirator.com



This product is manufactured under license to the 'Braun' trademark.
'Braun' is a registered trademark of Braun GmbH, Kronberg, Germany.

© 2018, All rights reserved. BNA 100US 25OCT18 P/N: A002262R1

Brand: Braun
Category: Nasal Aspirator
Model: **BNA100US**
Artwork Part #: A002262R1
Die Line Part #: NA
Subject: Quick Start Guide
Region: US
Flat Size: **MM: W375 x H90**
Folded Size: **MM: W75 x H90**
Scale: 1/1
Material: 160 gr/m2 Stock
Page count: 10
Revision: 4
Date: **25OCT18**
Release Date: **12NOV18**
Rerelease Date:

Colors:

Dielines (Do not print)

Cyan
100%

Magenta
100%

Yellow
100%

Black
100%

Spot Colors

PMS
645C



Coating:
Overall Gloss



Varnish:
Spot UV

Special Instructions:

Quality Requirement of Artwork and Quality
Clarification Process of Artwork Printing:
Meet Eng-QS-06&02

Helen of Troy
Creative Services
Marlborough, MA 01752 USA
+1 508 490 7000