

DO NOT PRINT LEGEND

**Helen
of Troy**

**Creative Department
Artwork Specifications**

Braun

Thermometers

BNT050EE

Artwork Part #: **A007753R2**

Dieline Part #: **N/A**

Subject: **Quick Start Guide**

Region: **EMEA**

Flat Size: **MM: W266.7 x H101.6**

IN: W10.5 x H4

Folded Size: **MM: W38.1 x H101.6**

IN: W1.5 x H4

Scale: **1/1**

Material: **120g SBS (art paper)**

of Pages: **1 page front and back**

Revision: **4**

Date: **30JUL24**

Release Date: **12AUG24**

Re-Release Date: **00MON00**

Colors

Cyan
100%

Magenta
100%

Yellow
100%

Black (K)
100%



Dielines do not print


Coating: **Coating**

Special Instructions:


Quality Requirement of Artwork
**Clarification Process of Artwork Printing: Meet
Eng-QS-06&02**

Helen of Troy Creative Services
Marlborough, MA 01752 USA
Questions? Contact your representative


1 Wykonywanie pomiaru




Należy celować w obszar między brwiami.



Maksymalna odległość.




Głowę należy trzymać nieruchomo.




Termometr należy utrzymywać w prawidłowej pozycji.



2 Wskazówki dotyczące pomiarów temperatury



Przed pomiarem: termometr przechowywać w temperaturze pokojowej przez 10 minut, a osoba poddawana pomiarowi powinna przebywać w pomieszczeniu przez 30 minut.



BRAUN


Sensian™ 3

Let's get started!


Zaczynamy!

BNT050


1 Taking a measurement




Aim between eyebrows.



Maximum distance.




Keep head still.




Keep thermometer in correct position.



2 Temperature taking tips



Before measurement: Keep thermometer in room temperature for 10 minutes and individual indoors for 30 minutes.



© 2025, All rights reserved.

Kaz Europe Sàrl
A Helen of Troy Company
Q-Center, Chaux 4
CH-1030 Bussigny – Switzerland

Certain trademarks used under license from The Procter & Gamble Company or its affiliates.

Sensian™ is a trademark of Helen of Troy Limited.

www.BraunHealthcare.com

Made in China.



REF BNT050EE



P/N: A007753R2
30JUL24